

Lab. 2

Microtechniques

“Histotechniques”






Sum of steps used for preparation of a
microscopic slide ready to examination

Harmand Ali

→ harmand.ali@ishiki.edu.iq

➤ What is histology?

Histology is the study of the microscopic anatomy (microanatomy) of cells and tissues of animals and plants.



Types of specimen

1. Autopsy

post-mortem specimens

2. Biopsy

specimens took from the living animals or humans...



Methods of preparing microscopic slides:


1. Sectioning method

Done by *Microtome*



2. Non-Sectioning method

Smear, squash, spreading



Processing for slide preparation

To perform histological examination, must first prepare tissues:

-  1 Fixation
-  2 Dehydration
-  3 Embedding
-  4 Sectioning/mounting
-  5 Staining
-  6 Microscopic examination



1- Fixation

- **Fixation** is a complex series of chemical events that comes by reaction between the **fixative** and **protein** which form a gel like subs.
- **The aim of fixation:**
 - 1- To prevent autolysis and bacterial attack.
 - 2- To prevent changing of the tissue's volume and shape during processing.
 - 3- prepare tissue for **clear staining** of sections.
 - 4- To get tissue as **close** as their **living state**.

Types of fixatives

- Formalin 10% (very common)
- Bouin fixative (expensive)
- Alcohols (70%)
- Gluteraldehyde (E.M)





Washing

- The process to get rid of the excessive fixative solution.
- Depend on the type of the fixative used the different washer are used.



2- Dehydration

- **Remove water from fixed tissue**
 - graded alcohol series followed by "clearing agent" (xylenes, histoclear, toluene) to remove alcohol (making tissue somewhat transparent).
 - **miscible** = liquids will mix to form a homogenous solution

Dehydration



50% EtOH / 50% Phosphate-buffered Saline (PBS)

75% EtOH / 25% PBS

90% EtOH / 10% PBS

95% EtOH / 5% PBS

100% EtOH

clearing agent

paraffin wax (if doing paraffin embedding/sectioning)

3- Clearing

- Replacing the dehydrating fluid with a fluid that is totally miscible with both the dehydrating fluid and the embedding medium.
- **miscible** = liquids will mix to form a homogenous solution

Some clearing agents:

Xylene

Toluene

Chloroform

Benzene

Choice of a clearing agent depends on:

1. Type of tissues to be processed
2. Safety factors.
3. Cost and convenience.
4. Speedy removal of dehydrating agent .
5. Minimal tissue damage .

4- Embedding




A process by which tissues are surrounded by a medium such as agar, gelatin, or wax which when solidified will provide sufficient external support during sectioning.

- **Why so ?**

- ❑ Structural support to tissue during sectioning
- ❑ Makes tissue easier to cut
- ❑ Tissue needs to withstand sectioning process
- ❑ Components must stay in natural positions

Embedding

Tissue embedded in "hard" Paraffin wax medium





5- Sectioning


- Using microtome to cut transparent thin slices from embedded tissue specimen & mount on microscope slides
- **Why slices should be very thin?**
 - Allows you to see internal structure of tissue.
 - Allows stains or specific markers such as antibodies to more easily infiltrate tissue.
 - Allows light to pass through tissue making structure visible.

Sectioning by microtome


- A **microtome** is a mechanical instrument used to cut biological specimens into very thin segments for microscopic examination.

Types:

- Rotary microtome
- Freezing microtome
- Sliding microtome
- Ultra microtome(E.M)




Copyright ©2006 by The McGraw-Hill Companies, Inc.
All rights reserved.




Sectioning by microtome

- We use microtome to cut transparent thin slices from embedded tissue specimens
 - generally **5-15** mm sections
 - Thickness depends upon tissue, subsequent procedures treatment.



Mounting sections

- A. 40° C water bath
 1. Flattens paraffin section
 2. Permits mounting on slide
- B. Gelatin & albumin
- C. Glass slides
- D. Oven / air dry





6- Staining

- Use many techniques to stain tissue sections w/ dyes
 - first must remove embedding medium, **rehydrate**(series **alcohol**) tissues (dyes generally in aqueous solution)
 - common stains - hemotoxylin & eosin (H&E), trichrome
- **Why staining is important?**
 - Creates higher contrast that allows observation of structures that are not visible in unstained tissues.
 - May reveal differences in chemical nature of certain regions of tissue.

Coverslipping / Mounting

mov



*A process in which we use a medium for adhesion of coverslip to the slide to protect the tissue from the dust, MO. & dirt.

Processs:

1. Remove the slides from the staining agent and dry it .
2. Coverslip using appropriate sized cover glass.
3. Label slides if necessary and arrange accordingly.



The slides are ready to be examined Microscopically



Thank you