Ministry of Higher Education and Scientific Research University of Ishik College of education Department of Biology

#### Practical Animal physiology Lab.4 (4th Grade)

#### By: Gashen Bakhtyar

Gashen.ahmed@gmail.com

## Packed Cell Volume (PCV) and Red Blood Cell Indices.

### **Objective:**

At the end of these laboratory students should be able to do PCV test, understand its procedure and techniques in the laboratory and hospitals.

- Whole blood is comprised of erythrocytes, platelets and leukocytes.
- ▶ The cells are suspended in the aqueous medium (plasma).
- In blood from healthy individuals, erythrocytes constitute the vast majority of cells.
- Plasma mainly consists of water (approximately 93%) but also of salts, various proteins and lipids as well as other constituents, e.g. glucose.

Hematocrit and PCV: Is defined as the proportion of red blood cells to plasma within a sample of blood. Following collection of the sample, the specimen is centrifuged. Due to their weight, the red blood cells are forced to the bottom of the test tube. A determination of the percentage of these packed cells in comparison to the plasma is then made.

#### Normal hematocrit values:

- At birth (newborn): 53-65%
- Six to 12 months: 33-40%
- Adult males: 42-52%
- Adult females: 35-47%
- Elderly: Slightly decreased
- Pregnancy: Decreased (dilutional)



#### Procedure

- 1. Lance your fingertip to obtain a drop of blood. Wipe away first drop of blood because it dilutes the sample with tissue fluid.
- 2. Fill the capillary tube (heparinized) two-thirds to three-quarters full with fingertip blood.
- 3. Seal one end of the tube with clay.
- 4. Place the filled tube in the micro-hematocrit centrifuge with the plug end to the outside.
- 5. Centrifuge at high speed (10000 rpm) for 5 minutes.
- 6. Read the hematocrit by placing the tube in the micro-hematocrit reader.

# Thank you